


User Manual for ACRoBaT device

designed by David Bjanès

dbjanès@uw.edu
copyright 2017



File List:

Bill of Materials: List of needed materials

Instruction Manual: Steps to assemble device, with listed part numbers from BOM

ACRoBaT: Published paper

Folder List:

code/arduino : c code to download on Arduino Pro Micro

code/logging-program: both source files and *.exe file for launching desktop logging program

build-files/pcb-brd-files: eagle schematic and board files for generating gerber files and/or ordering PCBs

build-files/stl-autocad-files: *.stl files for ordering 3D printed items from shapeways

videos/: several movies for animals during the task and various points throughout protocol

Notes:

If user wants to change LED cue functionality, modifying code below is necessary. Each Trial Type is found in file *ArduinoTrainerMicroPro0d_TrialTypes.ino*.

Change variables:

GBL_LED_Cues_Overt = true; // Controls whether LED is illuminated at the beginning of trial
GBL_LED_Cues_Hidden = false; // Controls whether LED is illuminated when target is entered

Below is a mapping of step numbers in the paper to actual values on the program (from file *ArduinoTrainerMicroPro0.ino*).

```
// *****  
// Trial Types  
// *****  
  
// Phase I // ACRoBaT Protocol Step  
const int TRIAL_SIPPER_20 = 1; // #1  
const int TRIAL_SIPPER_COUNTS_10 = 2; // #2  
const int TRIAL_TOUCH_JYSTK = 3; // #3  
  
// Phase II  
const int TRIAL_MOVE_JYSTK = 4; // #4  
const int TRIAL_MOVE_JYSTK_LEFT_NoTimeout = 5; // #5  
const int TRIAL_MOVE_JYSTK_RIGHT_NoTimeout = 6; // #5  
const int TRIAL_MOVE_JYSTK_CNTR_NoTimeout = 7; // #7  
const int TRIAL_MOVE_JYSTK_TOP_L_NoTimeout = 8; // (optional #7)  
const int TRIAL_MOVE_JYSTK_TOP_R_NoTimeout = 9; // (optional #7)  
  
// Phase III  
const int TRIAL_FIND_TARGET_OVERT_NoTimeout = 10; // #6 #8  
const int TRIAL_FIND_TARGET_OVERT_Timeout = 11;  
  
// Phase IV  
const int TRIAL_FIND_TARGET_HIDDEN_NoTimeout = 12; // #9  
const int TRIAL_FIND_TARGET_HIDDEN_Timeout_Step1 = 13; // #10  
const int TRIAL_FIND_TARGET_HIDDEN_Timeout_Step2 = 14; // #11  
const int TRIAL_FIND_TARGET_HIDDEN_Timeout_Step3 = 15; // #12  
const int TRIAL_FIND_TARGET_HIDDEN_Timeout_Step4 = 16; // #13  
const int TRIAL_FIND_TARGET_HIDDEN_Timeout_Step5 = 17; // #14  
const int TRIAL_FIND_TARGET_HIDDEN_Timeout_Step6 = 18; // #15  
const int TRIAL_FIND_TARGET_HIDDEN_EqualTimeout = 19; // #16
```